


3M Science.
Applied to Life.™



3M™ Impregum™ Super Quick
Polyether Impression Material

**Capture
every detail
in 2 minutes?
Yes, it can.**



Be impressed.

The fastest way to polyether precision for smaller cases

3M™ Impregum™ Super


Impregum Polyether impression material has always been the trusted solution for complex situations and implants. So wouldn't it be great if you could get that same outstanding impression accuracy and detail reproduction even in single-unit cases? To ensure that the highest demands are fulfilled without limiting productivity, we created new Impregum Super Quick Polyether impression material.

With its brand new formulation Impregum Super Quick material combines the outstanding performance and reliability of a polyether with the speed of a VPS. Deep margins and moist environments will no longer slow you down.



Quicker set.

So fast, so good: With just a two minute setting time the new Impregum Super Quick material allows for faster impressions in less complex cases like single-unit crowns, implants, inlays/onlays and smaller bridges up to 3 units.



100% Polyether. Trusted precision.

Like all polyethers Impregum Super Quick Materials demonstrate excellent performance in moist situations and great flow behaviour. Both lead to a high level of accuracy and detail reproduction for high-precision restorations.



It is especially important that you get the full working time of your impression material. Because it is a polyether, you get the full 45 seconds of working time with Impregum Super Quick material. You are able to avoid inaccuracies associated with premature setting that is experienced with VPS materials.



Make your patients feel comfortable.

An improved fresh, minty taste together with a short time in the mouth significantly increases your patients' comfort.



“I would, without a doubt, recommend it – once you've used Impregum, you won't use anything else.”

Dr Kunal Shah BDS
Principal of LeoDental, London


ether
es.
Quick.



Fast and gentle moisture management and retraction of gingival tissue with 3M™ Astringent Retraction Paste.



Convenient and highly accurate application of 3M™ Impregum™ Super Quick Light Body Polyether Impression Material using 3M™ Intra-oral Syringe Purple.



Precise, void-free details plus contrasting colors for excellent readability.

High-precision 1-step impression made with 3M™ Impregum™ Super Quick Polyether Impression Material (Medium and Light Body) and single-use 3M™ Impression Tray.

HIGHLY SATISFIED

A clinical evaluation shows a **very high satisfaction rate** of almost **95%** regarding the **precision** of 3M™ Impregum™ Super Quick Material impressions and the **fit of the final restorations.**



Source: 3M field evaluation with 447 participants from Europe and the US, 2017.

3M™ Impregum™ Super Quick

- ▶ **First superfast polyether with 45 sec. working time and quick 2-minute setting**
- ▶ **Available in medium and light body consistencies**
- ▶ **For monophasic and 1-step impression techniques**
- ▶ **Perfect for impressions for single-unit crowns, implants, inlays/onlays and smaller bridges up to 3 units**





	Dispensing system	Impression technique	Single-unit crown	Small implant (≤ 2 units)	Inlay/onlay (≤ 2 units)	Small bridge (≤ 3 units)	Recommended max. working time (min:sec)	Intra-oral setting time at 37°C (min:sec)
3M™ Impregum™ Penta™ Super Quick (Medium Body)		Monophase or 1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Penta™ Super Quick (Heavy Body)		1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Super Quick (Light Body)		1-step	✓	✓	✓	✓	0:45	2:00

Ordering information



3M™ Impregum™ Penta™ Super Quick Medium Body Introductory Kit (69383)



3M™ Impregum™ Penta™ Super Quick Medium Body/Light Body Introductory Kit (69381)



3M™ Impregum™ Penta™ Super Quick Medium Body Upgrade Kit (69404)



3M™ Intra-oral polyether syringes (71507)

Item no.	Product information
New Introductory kits for 3M™ Pentamix™ 3 Automatic Mixing Unit	
69383	3M™ Impregum™ Penta™ Super Quick Medium Body 1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 1 3M™ Penta™ Cartridge MB, 10 3M™ Penta™ Mixing Tips – Red, 1 3M™ Penta™ Elastomer Syringe, 1 Polyether Tray Adhesive – 17 ml, 3 3M™ Astringent Retraction Paste capsules
69381	3M™ Impregum™ Penta™ Super Quick Medium Body/Light Body 1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 2 cartridges LB – 50 ml each, 1 3M™ Penta™ Cartridge MB, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes Purple, 5 3M™ Garant™ Mixing Tips purple, 5 3M™ Garant™ Intra-oral Tips white, 1 Polyether Tray Adhesive – 17 ml, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Impression Trays
69413	3M™ Impregum™ Penta™ Super Quick Heavy Body/Light Body 1 base paste HB – 300 ml, 1 catalyst HB – 60 ml, 2 3M™ Impregum™ Super Quick LB cartridges – 50 ml each, 1 3M™ Impregum™ Penta™ HB Cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes Purple, 5 3M™ Garant™ Mixing Tips purple, 5 3M™ Garant™ Intra-oral Tips white, 1 3M™ Polyether Tray Adhesive – 17 ml, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Impression Trays
New Upgrade kits for 3M™ Pentamix™ 3 Automatic Mixing Unit	
69404	3M™ Impregum™ Penta™ Super Quick Medium Body 1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 10 3M™ Penta™ Mixing Tips – Red
69411	3M™ Impregum™ Penta™ 3 Super Quick Heavy Body/Light Body 1 base paste HB – 300 ml, 1 catalyst HB – 60 ml, 1 3M™ Impregum™ Super Quick LB cartridge – 50 ml, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple, 5 3M™ Garant™ Intra-oral Tips – White
New Refills	
69385	3M™ Impregum™ Penta™ Super Quick Medium Body Refill 2 base pastes - 300 ml each, 2 catalysts - 60 ml each
69406	3M™ Impregum™ Penta™ 3 Super Quick Heavy Body Refill 2 base pastes – 300 ml each, 2 catalysts – 60 ml each
69379	3M™ Impregum™ Super Quick Light Body Refill 4 cartridges – 50 ml each, 5 3M™ Intra-oral Syringes Purple, 5 3M™ Garant™ Mixing Tips purple
New Mixing devices and accessories	
77873	3M™ Pentamix™ 3 Automatic Mixing Unit UK
77902	3M™ Pentamix™ Lite Automatic Mixing Unit UK
69387	3M™ Impregum™ Penta™ Super Quick Cartridge for Pentamix™ 3
71507	3M™ Intra-oral polyether syringes (purple), 20 pieces
77580	3M™ Garant™ cartridge dispenser 1:1/2:1
56944	3M™ Astringent Retraction Paste capsules (25pk)
69413	Polyether Tray Adhesive – 17 ml
77949	3M™ Penta™ Mixing Tips – Red, 50 pieces
71454	3M™ Garant™ Mixing Tips purple, 50 pieces
71461N	3M™ Garant™ Intra-oral Tips white, 50 pieces
71616	3M™ Impression Trays, size M, upper jaw (10)
71619	3M™ Impression Trays, size M, lower jaw (10)

3M.co.uk/ImpSQ

3M Oral Care
3M United Kingdom PLC
Charnwood Campus
10 Bakewell Road
Loughborough
LE11 5RB

3M Ireland Ltd
The Iveagh Building
The Park, Carrickmines
Dublin D18 X015
Ireland
+353 (0)1280 3555

0800 626 578

OralCareUK@mmm.com

@3MOralCareUK

GD028. J442717.

3M, Garant, Impregum, Penta and Pentamix are trademarks of the 3M company. All other trademarks are owned by other companies. © 3M 2018. All rights reserved.

