

3M Science.
Applied to Life.™

3M™ RelyX™ Unicem and Unicem 2
Self-Adhesive Resin Cements

All-round reliability.



Cements:
Self Adhesive

Long-term
performer:
Cement



**3M™ RelyX™
Unicem**
Cements

Clinically proven self-adhesive resin cement


For simplicity, convenience and reliability, 3M™ RelyX™ Unicem Self-Adhesive Resin Cement is your 'go-to' product for a wide range of indications. Thanks to the dual-cure, self-adhesive chemistry, Unicem cements deliver high adhesion, resistance to discolouration, long-term stability and more.

Ideal for:

- ▶ general cementation
- ▶ oxide ceramic restorations
- ▶ posts and bridges

Results of *The Dental Advisor* 15 year clinical performance study

Results of 15-year recall of restorations cemented
with 3M™ RelyX™ Unicem/Unicem 2



Source: *The Dental Advisor* 15-year Clinical Performance Study, March 2017
Rating system: Excellent +++++

Formulated to meet the needs of your practice and your patients!

- ▶ Saves time by eliminating etching, priming and bonding steps
- ▶ Virtually no post-op sensitivity
- ▶ Strong
- ▶ Moisture tolerant
- ▶ Endo tips make cementing posts easy and reliable



3M™ Aplicap™ Capsule

- ▶ Ideally suited for inlays, onlays, crowns and posts
- ▶ Hygienic unidose
- ▶ Consistent mix with any mixing machine
- ▶ 3M™ RelyX™ Unicem Aplicap™ Elongation Tip for virtually void-free cementation of posts
- ▶ Working time: 2 minutes



3M™ Maxicap™ Capsule

- ▶ Ideally suited for multiple-unit bridges
- ▶ Hygienic unidose
- ▶ Consistent mix with any mixing machine
- ▶ Working time: 2.5 minutes



3M™ Automix Dispenser

- ▶ Suitable for all sizes of restorations
- ▶ Choose amount dispensed – 8.5g syringe gives approximately 15 applications
- ▶ Simple and time-saving one-step solution with no need for mixer, activators or appliers
- ▶ Comes with a regular tips for dispensing into crown, and wide tips



3M™ Unicem 2 Clicker™ Dispenser

- ▶ For all sizes of restorations
- ▶ Choose amount dispensed – 11g dispensed gives 40 applications
- ▶ Time-saving, easy and hygienic handling
- ▶ Delivers pre-measured doses for consistent and economical mixing
- ▶ No need for mixer, activators, appliers, mixing tips or devices

Recommended applications (++) Highly recommended / + Recommended)

Indication	Metal/ Metal based	Glass ceramics (including Li. Disilicate e.g. e.max®, VITA Mark II)	Oxide ceramics (e.g. 3M™ Lava™ Plus, Brux Zir®, Procera®)	Resin composites (e.g. VITAENAMIC®)
Inlays/Onlays	++	+	++	+
Crowns/Bridges	++	++	++	++*
Endodontic posts	++		++	++
Maryland bridges	+		+	
On implant abutments	++	+	++	++

*crowns only

Ordering Information

Item no. Product information

3M™ RelyX™ Unicem 2 Automix Refills		
56848	Shade Translucent	1 × 8.5g syringes
56846	Shade A2 Universal	10 x regular mixing tips
56847	Shade A3 Opaque	5 x wide mixing tips
		5 x endo tips
3M™ RelyX™ Unicem 2 Automix Mix Tips		
56851	30 mixing tips (regular)	
56918	15 mixing tips (wide); 15 intra-oral tips	
56919	15 mixing tips (wide); 15 endo tips	
3M™ RelyX™ Unicem 2 Clicker™ Dispenser Refills		
56874	Shade Translucent	1 × 11g Clicker dispenser
56875	Shade A2 Universal	1 x mix pad (70144)
56876	Shade A3 Opaque	Technique card
		Instructions for use
3M™ RelyX™ Unicem – Aplicap™ and Maxicap™ capsules		
56828	Trial pack 20 Aplicap capsules (10 x A2 Universal, 6 x Translucent, 4 x A3 Opaque)	
56829	Trial kit 20 Aplicap capsules (10 x A2 Universal, 6 x Translucent, 4 x A3 Opaque), Aplicap activator, Aplicap applier	
56814	Intro kit 50 Aplicap capsules (24 x A2 Universal, 16 x Translucent, 10 x A3 Opaque), Aplicap activator, Aplicap applier	

Item no. Product information

3M™ RelyX™ Unicem refills	
56815	50 Aplicap capsules (24 x A2 Universal, 16 x Translucent, 10 x A3 Opaque)
56816	50 Aplicap capsules, Translucent, 10 elongation tips
56817	50 Aplicap capsules, A3 Opaque, 10 elongation tips
56818	50 Aplicap capsules, A2 Universal, 10 elongation tips
56819	20 Aplicap capsules, A1
56820	20 Aplicap capsules, White opaque
3M™ RelyX™ Unicem – Maxicap™ capsules	
56830	Intro kit 20 Aplicap capsules (10 x A2 Universal, 6 x Translucent, 4 x A3 Opaque), Maxicap activator, Maxicap applier
3M™ RelyX™ Unicem – Maxicap™ capsules	
56831	20 Maxicap capsules (10 x A2 Universal, 6 x Translucent, 4 x A3 Opaque)
56832	20 Maxicap capsules, Translucent
56833	20 Maxicap capsules, A3 Opaque
56834	20 Maxicap capsules, A2 Universal

3M Oral Care
3M United Kingdom PLC
 Charnwood Campus
 10 Bakewell Road
 Loughborough
 Leicestershire LE11 5RB
 0845 602 5091

3M Ireland Limited
 The Iveagh Building
 The Park
 Carrickmines
 Dublin 18
 Ireland
 www.3M.co.uk/oralcare

GD025. J394376.
 e.max is a registered trademark of Ivoclar Vivadent.
 Brux Zir is a registered trademark of Glidewell Laboratories.
 Procera is a registered trademark of Nobel Biocare.
 Vita Enamic is a registered trademark of Vita Zahnfabrik.
 3M, Aplicap, Clicker, Lava, Maxicap
 and RelyX are trademarks of the 3M company.
 © 3M 2018. All rights reserved.

